

§ PicoIndenter 操作步驟

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

The figure shows the 'TriboScan Quasi' software interface. At the top, there are fields for 'Zero', 'Displ. (nm)', 'Actuation (μN)', and 'Sample (μN)', with values: Zero: z -523.8150, Displ. (nm): 18.6700, Actuation (μN): 95.0418. Below this is a 'Load Function' window with tabs for 'Quasi', 'Analysis', 'Calibration', and 'About'. The 'Quasi' tab is selected, and a red arrow points to it with the text '1. 切換到此標籤頁'. The main graph shows 'Displ (nm)' on the y-axis (ranging from -50.0 to 600.0) and 'Time (seconds)' on the x-axis (ranging from 0 to 20). The graph displays a red line that ramps up linearly from 0 to 600 nm at 10 seconds, then ramps down linearly back to 0 at 20 seconds. Below the graph is a 'Control Feedback' section with 'Normal Feedback Gains' (Proportional: 0.00000, Integral: 1.00000, Derivative: 0.00000, Adaptive: 0.00000) and checkboxes for 'Displ. Control (Idle Mode)', 'G-Control', and 'Video Graphics'. A red arrow points to the 'Tune Displ. PID' button with the text '2. 按下 Tune Displ. PID 鈕'. Other controls include 'Add Segment', 'Remove Segment', 'Peak Disp (nm): 600.00', and 'Perform Indent'.

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The screenshot displays the TriboScan Quasi software interface. At the top, there are status fields for Displ. (nm), Actuation (µN), and Sample (µN). The main window is titled 'Load Function' and contains a graph of Displ. (nm) vs. Time (seconds). The graph shows a red line increasing linearly from 0 to 600 nm over 10 seconds, followed by a black line decreasing linearly back to 0 nm over the next 10 seconds. Below the graph, there are several control panels:

- Segment # 1:** A table with parameters for the first segment, including Begin Disp (nm), End Disp (nm), Seg. Time (sec), Loading Rate (nm/s), Begin Time (sec), End Time (sec), # to Average, # of Segments, Data Points, Data Acq. Rate (pts/sec), and Total Number of Points.
- Peak Disp (nm):** A text box containing the value 600.00.
- Control Feedback:** A panel with 'Normal Feedback Gains' (Proportional, Integral, Derivative, Adaptive) and 'Feedback Mode' set to 'Displ. Control'. There are also checkboxes for 'Displ. Control (Idle Mode)', 'Control', 'G-Gain', and 'Video Graphics'.
- Buttons:** 'Perform Indent', 'Tune Displ. PID', and 'Graph To Clipboard'.

Four red arrows point to specific elements with blue text annotations:

1. 切換到此標籤頁 (Switch to this tab page) - points to the 'Quasi' tab.
2. 切換到 Displ. Control (Switch to Displ. Control) - points to the 'Displ. Control' radio button in the Feedback Mode section.
3. 設定好探針欲壓入樣品之行程及深度 (nm) (Set the travel and depth (nm) of the probe to be indented) - points to the 'Peak Disp (nm)' text box.
4. 按下 Perform Indent 鈕後會要求設定存檔路徑及檔名，即開始 indent (After clicking the Perform Indent button, you will be prompted to set the save path and filename, then start indent) - points to the 'Perform Indent' button.